





The 70th ASMS Conference, Minneapolis, MN

# Extremely Low Ion Counting Capabilities with Orbitrap FTMS

#### Konstantin O. Nagornov,<sup>1</sup> Natalia Gasilova,<sup>2</sup> Laure Menin,<sup>2</sup> Anton N. Kozhinov,<sup>1</sup> and <u>Yury O. Tsybin</u><sup>1</sup>

<sup>1</sup> Spectroswiss, Lausanne, Switzerland <sup>2</sup> Ecole Polytechnique Fédérale de Lausanne, Lausanne, Switzerland

Thursday, June 9, 2022: ThOE am Fundamentals: Unconventional Approaches in MS (Honoring R. Graham Cooks)

#### Ion Detection in Orbitrap FTMS

- Ion identity (m/z) is encoded as a frequency of ion oscillations in an ion trap
- Frequencies of ion oscillations are measured as time-domain signals (transients)
- Fourier transform (FT) decodes transients to reveal frequencies (m/z) values



#### Properties of Time-Domain Transients

- $R \sim T$  Resolution increases linearly with transient length, T
  - Longer transients provide higher resolution
- $S/N \sim z$  Sensitivity increases linearly with ion charge state, z
  - Multiply charged ions provide more abundant peaks in mass spectra

different of the first stand in the set of the stand over the stand over the set of the

Amplitude,

## $S/N \sim \sqrt{T}$ • Sensitivity increases as a square root of transient length, *T*

• Longer transients provide more abundant peaks in mass spectra

## $S/N \sim \sqrt{N}$ • Sensitivity increases as a sqrt of a number *N* of averaged transients

• Scan averaging (microscans) provide better quality mass spectra

#### Time-domain transients represent the truly unreduced (full) data

#### **Consequences of Time-Domain Transients Properties**

- Charge detection mass spectrometry (CDMS): single transient analysis
  - Individual multiply (highly) charged ions provide sufficiently strong ion signals
  - Increased length transients allow to detect even low charge state single ions Exploris: 4+ peptide ions, SNR = 3.5 with 1 s transients Makarov et al., IJMS, 2021, 116607
  - That enables Individual Ion Counting, I2MS, and related CDMS methods (e.g. DMT) Williams, Makarov, Kelleher, Heck, Chait, ...
  - Many applications benefit from averaging (many) transients and/or spectra
    - Complex mixture analysis, e.g., petroleomics Rogers, Marshall, ...
    - Top-down mass spectrometry for protein analysis, including mAbs
    - Intact mass measurements of large proteins, viruses

#### Conventionally, 100 – 1 000 transients are averaged

Unconventionally, what can be achieved when more transients are averaged?!

## Low Ion Counting in MS Applications

Rare singly (or multiply) charged ions can enter a mass spectrometer once every 10-100 seconds as isolated species or in a matrix of other, vastly more abundant, ions.



Ship diesel engine analysis, TP172

- Aerosol chemistry research
- Ambient air quality monitoring
- Complex mixture analysis
- Isotopic ratio analysis of elements
- Breath analysis

- Orbitraps offer high and ultra-high resolution capabilities
- Certain measurements can be done over an extended time long data averaging

#### How low ion counting capabilities can Orbitraps offer?

## Fullerene $(C_{60})$ – The Reference Molecule



## Fullerene (C<sub>60</sub>) Low Ion Counting – Single Scan Data



Simulations: Q Exactive HF,  $C_{60}$ , Tacq = 512 ms, single scan

## Fullerene (C<sub>60</sub>) Analysis – Transient Summation



**Rare events** 

Simulations: Q Exactive HF,  $C_{60}$ , Tacq = 512 ms, multiple scans

## Low Ion Counting – Transient Summation



#### SNR into Charges per Second Calibration

Detection period, s	SNR in 1 scan for z=+20	MSNR (Makarov SNR) for z=+1
0.192	1.78	0.089
0.384	2.57	0.1285
0.768	3.7	0.185
1.536	5.3	0.265
3.072	7.6	0.38

Denisov & Makarov, JASMS, 2009, 1486–1495

Charges/s = 
$$\frac{SNR}{\sqrt{N_{scans}} \cdot MSNR} / T_{acq}$$

A necessary number of **N**<sub>scans</sub> averaged scans to reach a given SNR





### **Experimental Set-Up**

- Extended duration measurements performed on a Q Exactive HF Orbitrap FTMS
- Time-domain transients acquisition with an add-on DAQ system (FTMS Booster)





#### **Ion Statistics**

APCI, SIM = 50 m/z, R = 120k, T<sub>acq</sub> = 256 ms, AGC = 2e5, IT<sub>max</sub> = 1ms total experimental time: 440 min



#### Experiment: C<sub>60</sub> Analysis with a Q Exactive HF





#### Uranium Isotopic Ratio Analysis: Extreme Dynamic Range



#### Detection of Rare Events (Singly Charged Ions)



Enabling isotopic ratio analysis for extremely challenging elements

### Conclusions

# Time-domain transient analysis enables ion detection over extended periods with only one charge injected into Orbitrap per 10 (...100) seconds measurement time

To detect ion signals with SNR > 4 and 512 ms time-domain transients:

- 1 charge per scan will require ~100 scans
- 1 charge per 10 scans ~5 000 10 000 scans
- 1 charge per 100 scans up to 1 million scans (suggested *in-silico*, to demonstrate!)
- The induced current detection realized in Orbitraps demonstrate its exceptional power
- Extended experiments reveal the depth of information contained in the time-domain transients
- Absorption mode FT mass spectra (unreduced) are readily available and provide equal performance to time-domain transients, even for the extremely low ion counting applications
- Understanding of these capabilities may benefit particular applications, including elemental analysis (isotopes), clinical, e.g., breath analysis, and environmental monitoring (matrix effects?)

#### The power and full capabilities of the Orbitraps (and ICRs!) are yet to be realized!







Richard Knochenmuss RKResearch LLC



R. Graham Cooks